



Pipet Filler, Safety Bulb

Product #:	1465100
AED Price:	Contact Hach
Available	

Unique cone adapter allows universal pipet filling. Simply squeeze rubber bulb to charge for filling, place pipet into cone adapter and insert into sample. Release the bulb to fill. Remove bulb and cover pipet end with finger for dispensing. Easy and safe. Autoclavable.

Specifications

Quantity:	1 each
What's included?:	Includes autoclavable elastic wristband.

What's included?

Includes autoclavable elastic wristband.