



Photometric Iodine (MEBAK) sample cuvette test, 25 tests

Product #:

LCK240

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

For Iodine values > 0.2

Specifications

According to standard:	MEBAK II
Description:	analogous MEBAK-Method
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	> 0.2
Method:	MEBAK Method
Number of tests:	25
Parameter:	Photometric Iodine sample (PIS)
Platform:	LCK
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)