



Hydrogen peroxide pipette test 1 - 10 g/L H₂O₂

Product #:

LCW058

Hazardous

AED Price:

Contact Hach



For determination of Hydrogen peroxide by the Peroxomolybdate method.

Specifications

Measuring range:	1 - 10 g/L H ₂ O ₂
Method:	Peroxomolybdate
Number of tests:	40
Parameter:	Hydrogen peroxide
Special remarks:	additional LCW906 round cuvettes needed
Storage conditions:	2 °C - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)