



AT1000 Series Potentiometric Titrator with 2 Burettes and 2 Pumps - Model AT1222

Product #:

AED Price: Call for ship date AT1222.98

Contact Hach

Automatic titration without all the complications

Simplify complex titration with the Titralab AT1000 from Hach. Straight from the box the AT1000 eliminates operator interpretation and manual processes you have come to expect with existing manual titration. The automatic system gives you more reliable results without complicated analysis. To make setup and titration simple, the AT1000 utilizes application specific kits to make it quick and easy for anyone to set up and operate a test.

Reliable titration results

Eliminates operator interpretation and manual processes with automatic titration to quickly deliver accurate and repeatable results.

Pre-programmed titration methods detect end points and eliminate manual calculations to make results easier to achieve without advanced programming.

Simple setup and titration

Application-specific functions to eliminate complex titration setup and analysis.

Hach's unique application kits make it quick for anyone to set up and operate a test.

Specifications

Application:	10
Automatic reagent addition:	Integrated peristaltic pumps for reagent addition
Beaker volumes accepted:	System specific beakers: Polypropylene, 50 mL & 150 mL with clips
	Generic beakers supported: Up to 250 mL low form glass
Burette for Titrant:	2
Burette motor resolution:	20,000 steps with electronic µstepping technology (128 µsteps/step)
Burette volumes available:	2.5/5/10/25 mL, ISO8655-3
Calibration:	Electrode & titrant (up to 5 buffers in Auto Mode, fixed & user defined)
Calibration curves display:	Direct & derivative titration curves, electrode calibration curve
Calibration intervals/alerts/reminder:	Yes, user defined period
Certifications:	Safety IEC/EN 61010-1; EMC IEC/EN 61326-1
Data logger:	On USB flash drive, CSV format, Excel compatible
Data storage:	Date, Time, Operator IDs, Sample IDs
Display type:	5,7"; Graphic colour; VGA
Electrode:	INTELLICAL (Plug & Play" digital), analogue, photocolorimetric; 2 inputs available"

Embedded procedures: Tubing bubble purging, reagent priming, burette and pump replacement Category II Installation category: English, German, French, Italian, Spanish Languages user interface: Maintenance scheduler: Yes, user defined Material: Latene (PP), silicone, hardened glass, POM, stainless steel Measurement method: Potentiometric (zero & imposed current), amperometric, colorimetric Model: AT1222 **Operating conditions:** 15 - 35 °C (59 - 95 °F), 20 - 80 % RH, non-condensing Operating humidity: 20 - 80 % (non-condensing) **Operating Interface:** Soft keypad (silicone) Operating temperature range: 15 - 35 °C (59 - 95 °F) Output: Electrode ports (x2), USB (x2), Serial (x1), Ethernet (x1) Parameter: mV/pH, Conductivity, Temperature (°C or °F) Password Protection: Yes, user defined protected elements 2 Pollution degree: Power requirements (Hz): 50/60 Hz 100/240 V AC Power requirements (Voltage): 100 mL/min Pump flow rate: 2 Pumps: **Required Application Package/s:** AP0005.AT1222 pH, Alkalinity & Hardness in Water AP0011.AT1222 pH, Total Acidity & Chlorides in Food & Beverage AP0013.AT1222 pH, Total Acidity, Free & Total SO, in Wine **Resolution:** Conductivity: ±0.5 % of reading Integrated, magnetic stirring, beakers up to 250 mL Sample stand: -5 - 40 °C Storage conditions: Supported peripherals: Printer, PC software, balance, sample changer, keyboard, mouse, barcode reader, sampling pump, propeller stirrer Titration modes: Sample, blank, sample with blank, QC sample, QC sample with blank Potentiometric (zero & imposed current), amperometric, colorimetric Titration types: Warranty: 24 months Weight: 4 kg What's included?: 1 ea. AT1222-Titralab AT1000 Titrator 1 ea. LZE144- Power Supply 1 ea. LZE141- 4 position tube holder 1 ea. LZE185- Electrode Storage Tube 1 ea. Basic User Manual Sold Separately: Don't forget to order your Application Pack and Reagents and consider purchasing an AS1000 Sample Changer.

1 ea. AT1222-Titralab AT1000 Titrator 1 ea. LZE144- Power Supply 1 ea. LZE141- 4 position tube holder 1 ea. LZE185- Electrode Storage Tube 1 ea. Basic User ManualSold Separately: Don't forget to order your Application Pack and Reagents and consider purchasing an AS1000 Sample Changer.

Required Accessories

- Titration Application Pack Water: pH, Alkalinity & Hardness (ISE) (Item AP0005.AT1222)
- Titration Application Pack Food & Beverage: pH, Total Acidity & Chlorides in Food & Beverage (Item AP0011.AT1222)
- Titration Application Pack, Food & Beverage: pH, Total Acidity & SO₂ in wine (Item AP0013.AT1222)