



FerroZine Iron Reagent Solution, 0.009-1.400mg/L Fe

Product #:

230149

Hazardous

AED Price:
Available

Contact Hach



For the determination of soluble ferric (Fe^{3+}) and ferrous (Fe^{2+}) Iron at trace levels.

Specifications

Instrument:	All except DR/850, DR/820, PC II
Measuring range:	0.009 - 1.400 mg/L Fe
Method:	FerroZine
Method number:	8147
Number of tests:	500 - 1000
Parameter:	Iron
Platform:	Solution
Quantity:	500 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)