



## Electrode Cleaning solution, RENOVO.X extra strong, 250 mL

**Product #:**

**S16M002**

**Hazardous**

AED Price:  
Available

Contact Hach



Electrode cleaning solution to remove difficult protein deposits. Convenient for pH sensor cleaning after use in dirty water or biological active samples.

---

### Specifications

Description:	When the solutions contain protein, for example, stronger cleaning agents are required.
Formula:	Sodium Hypochlorite
Parameter:	Cleaning solution
Quantity:	250 mL
Solution Type:	Cleaning