



## FerroZine Iron Reagent Solution Pillows, 0.009-1.400mg/L Fe

**Product #:**

**230166**

AED Price:  
Available

Contact Hach

**Hazardous**



For the determination of soluble ferric ( $\text{Fe}^{3+}$ ) and ferrous ( $\text{Fe}^{2+}$ ) Iron at trace levels.

FerroZine Iron Reagent determines soluble ferric ( $\text{Fe}^{3+}$ ) and ferrous ( $\text{Fe}^{2+}$ ) Iron. Digestion is needed for strongly chelated Iron complexes (such as those with EDTA) and Total Iron. For trace level Iron analysis. 50/pack

---

### Specifications

Instrument:	All except DR/850, DR/820, PC II
Measuring range:	0.009 - 1.400 mg/L Fe
Method:	FerroZine
Method number:	8147
Number of tests:	50
Parameter:	Iron
Quantity:	50 /pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)