



FerroMo Iron Reagent 1 Powder Pillows

Product #:

2543768

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Iron in presence of Molybdenum compounds.

Specifications

Measuring range:	0.01 - 1.80 mg/L Fe
Method:	FerroMo
Method number:	8365
Number of tests:	25
Platform:	Powder Pillows
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)