



FerroVer Iron Reagent, Accuvac, 0.02-3.00 mg/L Fe

Product #:

2507025

Hazardous

AED Price: Available







Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For determination of total dissolved Iron, as well as complexed iron compounds, by the 1,10 Phenanthroline Method.

Specifications

EPA compliant: Yes

Measuring range: 0.02 - 3.00 mg/L Fe

Method: FerroVer
Method number: 8008

Number of tests: 25

Parameter: Iron

Platform: Accuvac

Quantity: /pk

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)