



Ferrous Iron Reagent, Accuvac, 0.02-3.00 mg/L Fe (II)

Product #: 2514025

Hazardous

OBSOLETE ITEM
This item is no longer available.



Suggested replacements:

- 103769

Accuvac ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For determination of Iron in the ferrous form (Fe^{2+}) by the Phenanthroline method.

Specifications

Instrument:	All except PC II
Measuring range:	0.02 - 3.00 mg/L Fe^{2+}
Method:	1,10 Phenanthroline
Method number:	8146
Number of tests:	25
Parameter:	Iron
Platform:	Accuvac
Quantity:	25 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)