



## Molybdenum 1 Reagent Powder Pillows, 20 mL, pk/100

<b>Product #:</b>	<b>2352449</b>
AED Price:	Contact Hach
Available	

Powder reagent in sealed package for the determination of Molybdenum by the ternary complex method (Hach method 8169), 0.02 - 3.00 mg/L Mo. For 20 mL sample.

---

### Specifications

Method:	Ternary Complex
Number of tests:	100
Platform:	Powder Pillows
Quantity:	100 /pk

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)