



Nickel Reagent Set, Heptoxime Method Nickel Reagent Set, 0.02-1.80mg/L Ni

Product #:

2243500

Hazardous

AED Price:

Contact Hach

No availability shown



For the determination of Nickel, Method: Chloroform extraction using Heptoxime Indicator.

Specifications

EPA compliant:	Yes
Instrument:	All except DR/800, PC II
Measuring range:	0.02 - 1.80 mg/L Ni
Method:	Heptoxime
Method number:	8037
Number of tests:	50
Parameter:	Nickel

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)