



Lead reagent set, 5-150 μ g/L Pb

Product #: 2375000

AED Price: Contact Hach Available







Acid soluble Lead in a potable water sample is first concentrated on a Fast Column Extractor. The Lead is then eluted from the Extractor and determined colorimetrically with an indicator.

Specifications

Instrument: All except DR/800

Measuring range: $5 - 150 \mu g/L Pb$

Method: LeadTrak

Method number: 8317

Number of tests: 20

Parameter: Lead

Shipping Method: Lead Reagent set

Required Accessories

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)