



## QAC Reagent Set, Powder Pillows, 0.2-5.0mg/L CTAB

**Product #:**

**2459200**

**Hazardous**

AED Price:  
Available

Contact Hach



Powder reagents in sealed packages for the determination of Quaternary Ammonium Compounds (Direct Binary Complex Method).

CTAB = Cetyltrimethyl Ammonium Bromide

---

### Specifications

|                  |                                |
|------------------|--------------------------------|
| Instrument:      | All except DR/800, PC II       |
| Measuring range: | 0.2 - 5.0 mg/L as CTAB         |
| Method:          | Direct Binary Complex          |
| Method number:   | 8337                           |
| Number of tests: | 100                            |
| Parameter:       | Ammonium compounds, quaternary |
| Platform:        | Powder Pillows                 |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)