



Mercuric nitrate cartridge for digital titrator, 0.2256 N, 1 piece

Product #:

1439301

AED Price:
Available

Contact Hach

Hazardous



Ready-to-use titrant solution for chloride determination by Digital Titrator titration with mercuric nitrate, using diphenylcarbazone indicator. Hach method 8206. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifications

Measuring range:	10 - 160 mg/L Cl ⁻
Method:	Mercuric Nitrate
Method number:	8206
Number of tests:	40 - 100
Parameter:	Chloride
Platform:	Digital Titrator
Titrant:	0.2256 N Hg(NO ₃) ₂

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)