



Porphyrim 2 Reagent Powder Pillows

Product #:

2603649

Hazardous

AED Price:
Available

Contact Hach



Powder reagent in sealed package for the determination of Copper with Porphyrim method, 2-210 µg/l Cu.

For Copper determination by the Porphyrim method. Hach Method 8143 Range: 2 - 210 µg/L Cu. Sample Size: 10 mL Pack of 100 Powder Pillows.

Specifications

Method:	Porphyrim
Number of tests:	50 - 100
Platform:	Powder Pillows
Quantity:	100 /pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)