



Volatile Acids Reagent Set, 27-2800 mg/L HOAc

Product #:

2244700

Hazardous

AED Price:
Available

Contact Hach



For determination of Volatile Acids by the Esterification method.

Specifications

Instrument:	All except PC II
Measuring range:	27 - 2800 mg/L Acetic Acid
Method:	Esterification
Method number:	8196
Number of tests:	100
Parameter:	Volatile acids
Platform:	Solution
Special remarks:	Set contains: Ethylene glycol (203953), Ferric Chloride-Sulfuric Acid (204253), Hydroxyleamin Hydrochloride solution (81842), Sodium Hydroxide Standard solution 4.5 N (204053), Sulfuric Acid Standard solution 19.2 N (203832)
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)