



Tannin & Lignin Reagent Set, 0.1-9.0 mg/L

Product #:

2244600

Hazardous

AED Price:
Available

Contact Hach



For determination of Tannin and Lignin by the Tyrosine Method. For 25 mL sample.

Specifications

Instrument:	All except DR/820, PC II
Measuring range:	0.1 - 9.0 mg/L as Tannic Acid
Method:	Tyrosine
Method number:	8193
Number of tests:	100
Parameter:	Tannin & Lignin
Platform:	Solution

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)