



## Pipette validation kit

**Product #:** LCA722  
**AED Price:** Contact Hach

For accuracy check of pipettes within the Analytical Quality Assurance (AQA).

Errors in pipetted volumes are frequently attributable to a badly adjusted pipette. This error may occur after long use, overwinding, mechanical stress, or incorrect handling. Regular checks with the pipette test kit facilitate its early detection.

---

### Specifications

Description:	Pipette Test Kit
Method:	Decolourisation
Number of tests:	5
Parameter:	Pipette Test Kit
Platform:	LCA
Special remarks:	A provided colour solution of exactly known composition is diluted with the test volume from the Pipette tips. The resultant absorption is a measure of the deviation from the correct volume. Pipetting errors of more than 2% can be clearly detected in this way.
Storage conditions:	15 °C - 25 °C