



## DR6000 application software for brewery analysis

**Product #:**

**LZV936**

**AED Price:**

Contact Hach

The additional brewery analysis software is a compilation of all spectrophotometric applications that are relevant for brewery analysis.

The application software contains 12 brewery methods which can be transferred to the spectrophotometer via USB. The methods are based on published and observed brewery methods according to MEBAK and include procedures and calculations for: Anthocyanogens, Iron, Steam volatile phenols, Bitter units, Photometric iodine, Thiobarbituric acid number (TAN), Free amino nitrogen, Reductones, Vicinal diketones.