



Conditioning Solution, ISENH3181, 500 mL

Product #:	2541249
AED Price:	Contact Hach
Available	

Used as a conditioning solution between samples

Conditioning solution for Ammonia Ion Selective Electrode (ISENH3181), 500 mL bottle is used as a conditioning solution between samples.

This solution is for use between sample tests. It will help to reset the internal pH, and provide quicker, more accurate readings.

Please Note: This solution is **not** to be used for storage (overnight or longer), please see ISENH3181 electrode manual.

Specifications

Description:	ISENH3181
Formula:	$(\text{NH}_4)\text{HSO}_4$
Platform:	Solution
Quantity:	500 mL
Solution Type:	Storage