



## Singlet Single-Use pH Buffer, pH 10.01, pk/20

**Product #:** 2770220  
**AED Price:** Contact Hach  
**Available**

### Singlet Solution Packs

pH Buffer Standard Solution for pH electrode calibration. Sealed foil pouches of pH buffers. pH is  $\pm 0.02$  at 25°C. Pack of 20 of pH 10.01, 20 mL each

#### Ideal for on-site analyses

Sealed pouches are easy to transport - no risk of spillage.

#### Easy to Use

Just open the pouch and insert the electrode.

#### Fresh and reliable

Foil pouches prevent contamination and air exposure.

#### Design for use with Hach meters portfolio and other meters.

The basis for accurate pH measurements is good care and calibration of the pH Electrode.

---

### Specifications

|                     |                                 |
|---------------------|---------------------------------|
| Color coded?:       | Yes                             |
| Description:        | pH Buffer Calibration reagent   |
| NIST traceable:     | Yes                             |
| Package Type:       | Singlet                         |
| pH value(s):        | 10.01 $\pm 0.02$ pH at 25°C     |
| Quantity:           | 20 pcs                          |
| Special remarks:    | Discard solution after use      |
| Storage conditions: | 10 - 25 °C (protect from light) |