



Quantab Chloride Test Strips, Low Range, 30-600 mg/L, 40 pcs

Product #:

2744940

AED Price:
Available

Contact Hach

Ideal for determining quantitative results quickly, simply dip the strip in the sample and compare the length of the colour change to a printed scale.

Specifications

Measuring range:	30 - 600 mg/L Cl^-
Method:	Titration Strip
Number of tests:	40
Parameter:	Chloride
Platform:	Test Strip
Smallest Increments Steps:	10-20 ppm increments
Units:	40 tests