



## Filters, Hach 934-AH, 25 mm, 100/pk

<b>Product #:</b>	<b>2511203</b>
AED Price:	Contact Hach
Available	

Hach 934-AH Filters are specified for use in Standard Methods 2540 C, 2540 D and 2540 E. Fast flow rate, fine porosity and 1.5  $\mu\text{m}$  particle retention have made this binderless borosilicate glass microfiber filter the accepted industry standard for over 50 years. Hach 934-AH Filters can withstand temperatures over 500°C which is a unique and especially important feature for laboratories performing Fixed and Volatile Solids analysis of water and wastewater.

934-AH is a registered trademark of Whatman, Inc./GE Healthcare and used under license.

---

### Specifications

Diameter:	25 mm
Material:	Glass Microfiber
Pore size:	1.5 $\mu\text{m}$
Quantity:	100 /pk
Units:	100/pk