



Pipet, serological, sterile, polystyrene w/ plug, 10 mL, 500/pk

Product #: 2092868

OBSOLETE ITEM
This item is no longer available.

Disposable polystyrene with plug in mouth end. Calibrated "to deliver/blow out".

Sterile, disposable plastic pipet with plug. Bulk packaged (10x 50/pk). Color coded. Black ascending/descending graduations. Calibrated "To Deliver". Accuracy $\pm 2\%$.

Specifications

Capacity:	10 mL
Graduations:	0.10 mL
Quantity:	500 /pk