



## Thermometer, 30 to 50 C°, 250 mm

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2999200</b> |
| AED Price:        | Contact Hach   |
| Available         |                |

Non-mercury thermometer, 30 - 50 degrees C Length: 250 mm Immersion: 76 mm Safety coated.

Blue fill solution on white background provides excellent visibility.

Non-mercury thermometer, 30 - 50 degrees C Length: 250 mm Immersion: 76 mm Safety coated.

Blue fill solution on white background provides excellent visibility.

---

### Specifications

|                         |                   |
|-------------------------|-------------------|
| Length:                 | 250 mm            |
| Quantity:               | 1                 |
| Temperature range:      | 30 °C - 50 °C     |
| Thermometer divisions:  | 0.1 °C            |
| Thermometer fluid type: | Non-Mercury, blue |