



## Paddle tester, total aerobic bacteria, yeast and mold, 10 pieces

**Product #:**

**2610810**

AED Price:  
Available

Contact Hach

Hach Paddle Testers allow fast, easy screening of water samples, solid surfaces, and nonviscous liquids where unsanitary conditions may exist. The paddle is a double-sided slide attached to the vial cap. Each side of the slide is used to perform a separate test. The slide has a molded-in grid for easy colony counting without opening the optically clear, leakproof vial. Paddle Tester, Total Aerobic Bacteria/Yeast and Mold (Tryptone Glucose Extract Agar/Sabouraud Dextrose Agar), pk/10.

### **Save time and simplify lab and field testing with Hach prepared P/A media**

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### **Report Total Coliform values with EPA accepted Presence-Absence bottles and media**

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### **Simplify routine testing with self-contained, easy-to-use Paddle Testers and BART products**

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## Specifications

|             |                                        |
|-------------|----------------------------------------|
| Media Type: | Total Aerobic Bacteria, Yeast and Mold |
| Quantity:   | 10 /pk                                 |