



Nitrate Reagent Powder Pillows, 0.3-30mg/L NO₃-N

Product #:

2106169

AED Price:
Available

Contact Hach

Hazardous



Powder reagent in sealed packages for the determination of Nitrate by Cadmium reduction method. For 10 mL sample, 100 pieces.

Specifications

Description:	NitraVer 5 Nitrate Reagent Powder Pillows, 10mL HR, MR
Instrument:	All except DR/3000, DR/2010, DR/2000
Measuring range:	0.3 - 30.0 mg/L NO ₃ -N
Measuring range (2):	0.1 - 10.0 mg/L NO ₃ -N
Method:	Cadmium Reduction
Method number:	8039 HR
Method number :	8171 MR
Number of tests:	100
Parameter:	Nitrate
Platform:	Powder Pillows
Quantity:	100 /pk
Reference method:	SM 4500-NO3-B
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)