



Nitrate TNT, 0.2 - 30 mg/L NO₃-N

Product #:

2605345

AED Price:

Contact Hach

Ships within 2 weeks

Hazardous



Reagent Vials (Ø 16 mm) contain premeasured DPD reagent, optimised for reliable measuring results and easy handling. TNT = Test 'N Tube

Specifications

Description:	NitraVer X Nitrate Reagent Set
Instrument:	DR 900
Measuring range:	0.2 - 30.0 mg/L NO ₃ -N
Method:	Chromotropic Acid
Method number:	10020
Number of tests:	50
Parameter:	Nitrate
Platform:	Test 'N Tube
Quantity:	50 /pk
Special remarks:	Set contains: Cuvette test tubes, NitraVer X Reagent A (2605400), NitraVer X Reagent B Powder Pillows (2605546)
Storage conditions:	10 °C - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)