



PHG311-9 Glass pH Electrode for Alkaline Samples, with screw cap, without cable

Product #:	E11M004
AED Price:	Contact Hach
Available	

Glass pH electrode with screw cap for alkaline samples. Needs cable.

Specifications

Compatible instruments:	MeterLab & sensION+ GLP lab
Connector:	Screw cap
Electrode type:	PHG311-9
Length:	120 mm
Material:	Glass
Measuring range:	0 - 14 pH
Parameter:	pH
Sensor type :	Glass (for alkaline samples)
Special Features:	Glass membrane for reduced sodium error
Temperature range:	0 - 80 °C
Thermistor:	none
Type:	Lab, external reference probe needed
Warranty:	6 months
What's included?:	Includes: Radiometer Analytical PHG311-8 pH electrode with Certificate of Conformity.

What's included?

Includes: Radiometer Analytical PHG311-8 pH electrode with Certificate of Conformity.

Required Accessories

- CL114 cable, 1m, for screw cap electrode FX/S7/coax, BNC plug (Item A94L114)
- REF321 Reference electrode, Ag/AgCl (Item E21M002)