



PHC3001 Combined pH electrode, Ag/AgCl, BNC

Product #:	E16M300
AED Price:	Contact Hach
Available	

General-purpose pH electrode with Ag/AgCl reference system. Measuring range 0-12 pH, 0-80°C. Liquid junction: porous pin. BNC plug.

Specifications

Cable length:	1 m
Compatible instruments:	MeterLab & sensION+ GLP lab
Connector:	BNC Connector
Dimensions (D x L):	12 mm x 103 mm
Electrode junction :	Porous pin
Electrode type:	Laboratory combination pH electrode
Filling solution:	Refillable: 3 M KCl + sat. Ag/Cl
Length:	103 mm
Material:	Glass
Measuring range:	0 - 12 pH
Parameter:	pH
Sensor type :	Glass
Temperature range:	0 - 80 °C
Thermistor:	none (external required)
Type:	Lab, combined electrode

Required Accessories

- Filling Solution, Reference, 3 M KCl with AgCl, 100 mL (Item S21M004)