



Test Kit Iron, 0 - 4 mg/L Fe, Colour disc method, 100 tests

Product #:

146400

Hazardous

AED Price:
Available

Contact Hach



Specifications

Case style:	B
Measuring range:	Medium Single Range
Method:	Colour disc/ 1,10 Phenanthroline
Model:	IR-18
Number of tests:	100
Parameter:	Iron
Platform:	Color Disc
Smallest increment:	0.1
Standard range:	0 - 4 mg/L Fe
Test kit model:	IR-18
Type:	Single Metric