



Test Kit, chlorine (free & total) & pH, model CN-67

Product #:

1411100

Hazardous

AED Price:
Available

Contact Hach



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

One of Hach's most popular test kits

Includes color discs and a comparator for monitoring chlorine and pH

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Application:	Pool and Spa
Measuring range:	6.6 - 8.4 pH units
Measuring range (2):	0 - 3.4 mg/L
Model:	CN-67
Number of tests:	300
Parameter:	Chlorine (Free & Total) and pH
Platform:	Color Disc