



## Hardness, iron, and pH test kit, model HA-62A

**Product #:**

**183701**

AED Price:  
Available

Contact Hach

**Hazardous**



Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

**Popular tests at a great price in Hach's convenient and protective test kit carrying case.**

**Uses color disc and drop count titration methods**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Application:	Water Quality
Measuring range:	6.6 - 8.4 pH units
Model:	HA-62A
Number of tests:	100
Parameter:	Hardness, Iron, and pH (6.6-8.4)
Platform:	Color Disc