



TEST KIT, 5-EP mg/L HARDNESS

Product #:

145401

Hazardous

AED Price:
Available

Contact Hach



Total Hardness drop count titration/EDTA test kit with range of 20 - 400 mg/L; contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

| | |
|---------------------|--|
| Case style: | B |
| Measuring range: | 20 - 400 mg/L CaCO_3 |
| Method: | Drop count titration/ EDTA-powder pillows |
| Number of tests: | 100 |
| Parameter: | Hardness, total - As CaCO_3 |
| Platform: | Drop Count Titration |
| Smallest increment: | 20 |
| Test kit model: | 5-EP MG-L |