



Nitrogen, Nitrite Colour Disc Test Kit, Model NI-6

Product #:

224000

Hazardous

AED Price:
Available

Contact Hach



Utilising the simple colour disc method/ferrous sulfate, this kit measures Nitrogen, Nitrite as NO_2^- and contains 100 tests. Range (mg/L): 0-8/0-1600.

Simple, low-cost colour disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifications

Case style:	D
Measuring range:	Single Range
Method:	Colour disc/ Ferrous sulphate
Number of tests:	100
Parameter:	Nitrogen, Nitrite
Platform:	Color Disc
Smallest increment:	2 40
Test kit model:	NI-6
What's included?:	Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.

What's included?

Includes colour disc, comparator box and tubes, reagents, deionised dilution water, instructions and carrying case.