



Cylinder, graduated, mixing, glass, 50 mL ± 0.4 mL, 1.0 mL divisions, glass stopper #16

Product #:	189641
AED Price:	Contact Hach
Available	

Tall form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Specifications

Capacity:	50 mL
Certifications:	Class B
Housing material:	Glass
Measurement method:	To Contain
Quantity:	1 each
Type:	Single Metric, Mixing