



Beaker, PP, low form, 100 mL

Product #:	108042
AED Price:	Contact Hach
Ships within 3 days	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Capacity:	100 mL
Diameter:	63.0 mm with spout
Housing material:	Polypropylene
Material:	Polypropylene, Low Form
Quantity:	1 each
Units:	1 each