



Nitrification inhibitor for BOD (35 g)

Product #: 253335
AED Price: Contact Hach
Available

TCMP - 2-chloro-6 (trichloromethyl) pyridine (N-Serve), coated on an inert substrate. Use 0.16 gr in a 300-mL sample inhibits nitrogenous oxygen demand in the BOD test (Use Dispenser Cap, Part No. 45901), 35 g.

Specifications

Application:	Nitrification Inhibitor (CBOD)
Description:	Nitrification Inhibitor
Formula:	Trichloromethyl Pyridine (TCMP)
Method number:	8043
Number of tests:	200
Package Type:	Bottle
Parameter:	Biochemical Oxygen Demand
Platform:	Jar
Quantity:	35

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)
- Dispenser cap, formula 2533 nitrification inhibitor 35 g bottle, 1/pk (Item 45901)