



Total Chlorine Reagent Set for chlorine analyser CL17/CL17sc

Product #:

2557000

AED Price:
Available

Contact Hach

Hazardous



Formulated to meet the strictest quality standards

The Total Chlorine Reagent set is used with the Hach CL17/CL17sc analyzers as well as bulk analytical methods. Hach reagents provide a guarantee of accurate, consistent results.

Total Chlorine Reagent Set Consists of:

- One 24-gram bottle of DPD compound
- One 473-mL bottle of total chlorine indicator solution
- One 473-mL bottle of total chlorine buffer solution

*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

Specifications

| | |
|---------------------|---|
| Instrument: | CL17/CL17sc |
| Measuring range: | 0.03 - 5.00 mg/L Cl_2 for CL17, 0.03 - 10.00 mg/L Cl_2 for CL17sc |
| Method: | DPD |
| Method number: | 8370 |
| Number of tests: | 450 |
| Parameter: | Chlorine, Total |
| Storage conditions: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |
| What's included?: | DPD compound, 24 g |
| | Total chlorine indicator solution, 473 mL |
| | Total chlorine buffer solution, 473 mL |

What's included?

DPD compound, 24 gTotal chlorine indicator solution, 473 mLTotal chlorine buffer solution, 473 mL

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)