



BioTector B3500 Reagent, 6.0 N Sulfuric Acid with Mn Catalyst, 20 L

Product #:

2038162

Hazardous

AED Price:

Contact Hach



Specifications

Instrument: B3500

Measuring range: N/A

Parameter: TOC

Storage conditions: Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)