



KNF300 Pump 12 Month Service Kit

Product #: 19-KIT-111
AED Price: Contact Hach

For BioTector Series 4 analysers & B7000 supplied before July 2014 – service kit to be used every 12 months for the maintenance of the KNF 300 circulation pump

Specifications

What's included?: Includes: PTFE Diaphragm, NF300 O-Ring, Intermediate Plate, Connection Plate and PTFE Tape

What's included?

Includes: PTFE Diaphragm, NF300 O-Ring, Intermediate Plate, Connection Plate and PTFE Tape