



pH Electrode for 8362 Ultrapure Water Sensor

Product #: 08362=A=2000

AED Price: Contact Hach

Part of analogue pH sensor 8362, specially designed for pure water pH measurement.

Specifications

Application:	Ultrapure water
Cable connection:	S7 connector
Connector:	S7 Connector
Length:	120 mm
Material:	Sensor: Glass
Measuring range:	0 - 14 pH
Reference electrode:	None
Sensor thread:	Pg 13.5
Special notes:	pure water pH electrode
Temperature sensor:	PT100
Warranty:	6 months