



SOTA10 pH electrode, 10 m fixed cable, without T-sensor

Product #: LZX587
AED Price: Contact Hach

Standard glass and reference system

Specifications

Cable length:	10 m fixed cable
Diaphragm:	Ring diaphragm, porous PTFE, double
Electrolyte:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Immersion depth:	105 mm
Length:	140 mm
Measuring range:	0 - 14 pH
Operating temperature range:	0 - 80 °C
Pressure range:	3.5 bar
Warranty:	6 months