



Replacement ORP electrode for 1200-S sc sensor, with fitting

Product #: LZX890
AED Price: Contact Hach

Specifications

Accuracy:	± 0.2 °C
Cable connection:	S7 connector
Calibration method:	One point, standard buffer solution (automatic), One point, comparative measurement
Diameter:	12 mm
Flow:	max. 4 m/s
Length:	120 mm
Measuring range:	-1500 - 1500 mV
Operating temperature range:	-5 - 50 °C
Pressure range:	2 bar
Repeatability:	± 0.5 %
Response time:	< 2 min temperature
Storage conditions:	-20 °C - 60 °C
Warranty:	6 months
Weight:	< 0.15 kg