



## Iron (II / III) cuvette test 0.2-6.0 mg/L Fe, 24 tests

**Product #:**

**LCK320**

**AED Price:**

Contact Hach

**Hazardous**



### **LCK cuvette tests: Minimum amounts of reagents, maximum safety.**

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### **Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### **Specifications**

According to standard:	DIN 38406-E1
Description:	Iron II and III
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.2 - 6.0 mg/L Fe (II/III)
Method:	1.10-Phenanthroline
Number of tests:	24
Parameter:	Iron
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)

---

### **Required Accessories**

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)