



Cadmium cuvette test 0.02-0.3 mg/L Cd, 25 tests

Product #:

LCK308

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.02 - 0.3 mg/L Cd
Method:	Cadion
Number of tests:	25
Parameter:	Cadmium
Platform:	LCK
Special remarks:	For removal of interfering Calcium and/or Magnesium ions (up to 5000 mg/l) use LCW903 Calcium separation set before Cadmium determination with LCK308. For total Cadmium digestion use CRACK SET LCW902.
Storage conditions:	2 - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)