



Chromium (III and VI) cuvette test 0.03-1.0 mg/L Cr, 25 tests

Product #:

LCK313

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	EN ISO 11083, DIN 38405-D24
Description:	Chromium VI and Total
Digestion required:	Yes, for Total Chromium
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.03 - 1.0 mg/L Cr VI
Measuring range (2):	0.03 - 1.0 mg/L Cr total
Method:	Diphenylcarbazide
Number of tests:	25
Parameter:	Chromium
Platform:	LCK
Special remarks:	Test contains digestion for Total Chromium
Storage conditions:	2 - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)