



## Potassium cuvette test 8-50 mg/L K, 24 tests

**Product #:**

**LCK328**

**Hazardous**

**AED Price:**

Contact Hach



### **LCK cuvette tests: Minimum amounts of reagents, maximum safety.**

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

#### **Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### **Specifications**

EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	8 - 50 mg/L K
Method:	Kalignost
Number of tests:	24 (test includes zero solution)
Parameter:	Potassium
Platform:	LCK
Storage conditions:	15 - 25 °C

---

### **Required Accessories**

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)