



Silver cuvette test 5-2500 mg/L Ag, 25 tests

Product #:

LCK355

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

Description:	Two measuring ranges: (I) 5-400 mg/L Ag, (II) 400-2500 mg/L Ag
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	5 - 400 mg/L Ag (I)
Measuring range (2):	400 - 2500 mg/L Ag (II)
Method:	Hach Method
Number of tests:	25
Parameter:	Silver
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)